

Product Datasheet

Rat GHRL ELISA Kit (orb781589)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9QYH7
Sensitivity	28 pg/mL
Expiration Date	Please enquire.

