

Product Datasheet

Other HA ELISA Kit (orb781706)

| | |
|----------------------------|---|
| Description | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| Reactivity | All |
| Range | 1.57-100 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Sensitivity | 0.51 ng/mL |
| Expiration Date | Please enquire. |

