


## Product Datasheet

### Human Anti-HPF4 ELISA Kit (orb781874)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Human
<b>Range</b>	1.57-100 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates and other biological fluids.
<b>Sensitivity</b>	0.59 ng/mL
<b>Expiration Date</b>	Please enquire.



Standard curve of Human Anti-HPF4 ELISA

