

Product Datasheet

Sheep IL17 ELISA Kit (orb781937)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Sheep
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9TUP4
Sensitivity	5.6 pg/mL
Expiration Date	Please enquire.

