

Product Datasheet

Porcine IL8Ra ELISA Kit (orb782264)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Porcine |
| Range | 0.16-10 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | tissue homogenates, cell lysates and other biological fluids |
| Sensitivity | 0.054 ng/mL |
| Expiration Date | Please enquire. |

