

## Product Datasheet

### Human BTLA ELISA Kit (orb782463)

<b>Description</b>	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
<b>Reactivity</b>	Human
<b>Range</b>	0.47-30 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, urine, tissue homogenates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	BTLA
<b>Sensitivity</b>	0.22 ng/mL
<b>Expiration Date</b>	Please enquire.

