


Product Datasheet

Canine SDMA ELISA Kit (orb782628)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Canine
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Sensitivity	0.049 ng/mL
Expiration Date	Please enquire.



Standard Curve of Canine SDMA ELISA Kit.

