

## Product Datasheet

### Fish GSH ELISA Kit (orb782655)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Zebrafish
<b>Range</b>	1.57-100 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Sensitivity</b>	0.48 ng/mL
<b>Expiration Date</b>	Please enquire.

