

## Product Datasheet

### Human CP ELISA Kit (orb782665)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 0.16-10 ng/mL   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids    |
| <b>Sensitivity</b>         | 0.048 ng/mL   |
| <b>Expiration Date</b>     | Please enquire.   |

