

Product Datasheet

Rat p-hydroxyphenyllactic acid ELISA Kit (orb782666)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Rat
Range	15.6-1000 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Sensitivity	6.65 ng/mL
Expiration Date	Please enquire.

