

## Product Datasheet

### Rat TBARS ELISA Kit (orb782699)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| <b>Reactivity</b>          | Rat   |
| <b>Range</b>               | 6.25-400 nmol/mL  |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Sensitivity</b>         | 3.55 nmol/mL  |
| <b>Expiration Date</b>     | Please enquire.   |

