

## Product Datasheet

### Rat 5-HT ELISA Kit (orb782748)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Rat
<b>Range</b>	14.07-900 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Sensitivity</b>	3.97 ng/mL
<b>Expiration Date</b>	Please enquire.

