

## Product Datasheet

### Mouse Testo ELISA Kit (orb782823)

<b>Description</b>	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
<b>Reactivity</b>	Mouse
<b>Range</b>	62.5-4000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Sample Types</b>	serum, plasma, other biological fluids
<b>Sensitivity</b>	17.7 pg/mL
<b>Expiration Date</b>	Please enquire.

