


Product Datasheet

Human NPT(Neopterin) ELISA Kit (orb782859)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Human
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Sensitivity	0.048 ng/mL
Expiration Date	Please enquire.



Standard Curve of Human Neopterin ELISA ...

