

Product Datasheet

EPB41 (phospho-Tyr418 + Tyr660) antibody [Out of stock] (orb783184)

Description	EPB41 (phospho-Tyr418 + Tyr660) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human EPB41 around the phosphorylation site of Tyr418 + Tyr660 (NI(p-Y)IR)
Target	phospho-EPB41 (Tyr418 + Tyr660)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	97 kDa
Uniprot ID	P11171
Expiration Date	12 months from date of receipt.