

Product Datasheet

Mu Opioid Receptor (phospho-Ser375) antibody [Out of stock] (orb783186)

Description	Mu Opioid Receptor (phospho-Ser375) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from rat mu Opioid Receptor around the phosphorylation site of Ser375 (HP(p-S)TA)
Target	Phospho-mu Opioid Receptor (Ser375)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	P35372
Expiration Date	12 months from date of receipt.