
Product Datasheet

C5a Anaphylatoxin antibody (HRP) (orb783530)

Description

Rabbit polyclonal antibody against C5a

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse C5a (725-770aa)
Target	C5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	9/113kDa
Uniprot ID	P02748
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.