

Product Datasheet

MERTK antibody (orb783575)

Description	Mouse monoclonal antibody to MERTK
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human MERTK protein
Target	MERTK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	hIgG4(S228P)
Clonality	Monoclonal
Uniprot ID	Q12866
Expiration Date	12 months from date of receipt.