

Product Datasheet

Claudin-18.2 antibody (orb783578)

Description	Mouse monoclonal antibody to Claudin-18.2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Recombinant Human Claudin-18.2 protein
Target	CLDN18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	hIgG1
Clonality	Monoclonal
Uniprot ID	P56856
Dilution Range	Recommended dilution: FC:1:50-1:200
Expiration Date	12 months from date of receipt.