

Product Datasheet

Novel Coronavirus Nucleoprotein (N) antibody (orb783580)

Description	Mouse monoclonal antibody to Novel Coronavirus Nucleoprotein (N)
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, GICA
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa)
Target	N
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse scFv fusion with human IgG1 Fc
Clonality	Monoclonal
Uniprot ID	PODTC9
Dilution Range	ELISA:1:10000-1:50000, LF:1:500-1:5000
Expiration Date	12 months from date of receipt.