

## Product Datasheet

# Novel Coronavirus Spike glycoprotein(S) antibody (Biotin) (orb783581)

<b>Description</b>	Mouse monoclonal antibody against Novel Coronavirus Spike glycoprotein(S) conjugated to Biotin
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Human Novel Coronavirus Spike glycoprotein(S) (319-541aa)
<b>Target</b>	S
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	VHH fusion with human IgG1 Fc
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>PODTC2</b>
<b>Dilution Range</b>	Recommended dilution: ELISA:1:10000-1:50000
<b>Expiration Date</b>	12 months from date of receipt.