


## Product Datasheet

### Upf1/Rent1 (phospho-Thr28) antibody (orb7836)

<b>Description</b>	Rabbit monoclonal antibody to Upf1/Rent1 (phospho-Thr28)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IB, WB
<b>Immunogen</b>	Synthetic peptide linked to KLH
<b>Target</b>	Upf1/Rent1 (pThr28) [P-CS-1]
<b>Preservatives</b>	Purified antibody containing PBS + 0.09% sodium azide
<b>Form/Appearance</b>	Purified antibody containing PBS + 0.09% sodium azide
<b>Concentration</b>	0.16mg/ml
<b>Storage</b>	Store at +4 °C for 1 month, for long-term storage store at -20 °C
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Peptide affinity chromatography
<b>Uniprot ID</b>	<a href="#">Q92900</a>
<b>Entrez</b>	<a href="#">5976</a>
<b>Dilution Range</b>	Optimal antibody dilution should be determined by titration, however as a guideline start at 1µg/ml for western blot
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of using Upf1/Rent1...