


Product Datasheet

ICAM4 antibody (orb79346)

Description	Goat polyclonal antibody to ICAM4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ICAM4 (aa169-182)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 28kDa band observed in Human Umbilical Cord and in Erythrocytes lysates (calculated MW of 29.3kDa according to NP_001535.1). Recommended concentration: 1-3 µg/ml. Experiment Notes: This antibody is expected to recognize isoform 1 (NP_001535.1) and isoform 2 (NP_071772.1) only.
Protein Sequence	RHGSRVIYSESLER
Clonality	Polyclonal
MW	29.3; 25.9; 29.6
Entrez	3386
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Umbillica...