

Product Datasheet

ZHD11 antibody (orb799043)

| | |
|----------------------------|---|
| Description | ZHD11 antibody |
| Species/Host | Rabbit |
| Reactivity | A. thaliana |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Arabidopsis thaliana ZHD11 protein |
| Target | ZHD11 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9SEZ1 |
| Expiration Date | 12 months from date of receipt. |