


Product Datasheet

Human Epigen protein (orb80017)

Description	Epigen Recombinant Human produced in E.coli is a single, non-glycosylated, polypeptide chain contain...
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	EPGN
Preservatives	EPGN was lyophilized from 20 mM PBS buffer pH-7.4C ontaining 130mM sodium chloride.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	AVTVTPPITA QQADNIEGPI ALKFSLCLE DHNSYINGA CAFHHELEKA ICRCFTGYTG ERCEHLTLTS YA.
Purity	> 98.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized Epigen in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Epigen protei...