


Product Datasheet

Human IL33 protein (orb80083)

Description	Human IL33 Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain containin...
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	IL33
Preservatives	The protein is formulated in 20 mM Tris pH-8.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MSITGISPIT EYLASLSTYN DQSITFALED ESYEIVEDL KKDEKKDKVL LSYYESQHPS NESGDGVDGK MLMVTLSPK DFWLHANNKE HSVELHKCEK PLPDQAFFVL HNMHSNCVSF ECKTDPGVFI GVKDNHLALI KVDSSNLCT ENILFKLSET.
Purity	> 95.0% as determined by RP-HPLC and analysis by SDS-PAGE
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL33 protein