

Product Datasheet

Human MIF protein (orb80089)

Description	Human MIF His tag C-Terminus protein
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	Human MIF
Preservatives	Human MIF was lyophilized from 1mg/ml solution containing PBS pH-7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MPMFIVNTNV PRASVPDGFL SELTQQLAQA TGKPPQYIAV HVVPDQLMAF GGSSEPALCS LHSIGKIGGA QNRSYSKLLC GLLAERLRIS PDRVYINYD MNAANVGWNN STFALEHHHH HH.
Purity	> 95.0% as determined by(a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.
Solubility (25°C)	It is recommended to reconstitute the lyophilized MIF in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.