


Product Datasheet

Human TNF alpha His tag protein (orb80108)

Description	Human TNF alpha His tag protein
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	Human TNF alpha His tag
Preservatives	TNF-a His is supplied in 1xPBS, 50% glycerol.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95.0% as determined by RP-HPLC and analysis by SDS-PAGE
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human TNF alpha His...