


Product Datasheet

Human Thymic Stromal Lymphopoietin protein (orb80113)

Description	TSLP Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	TSLP
Preservatives	Lyophilized from concentrated solution (1mg/ml) contains 130mM NaCl, and 20 mM sodium phosphate, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	MYDFTNCDFE KIKAAYLSTI SKDLITYMSG TKSTEFNNTV SCSNRPHCLT EIQLTFNPT AGCASLAKEM FAMKTKAALA IWCPGYSETQ INATQAMKKR RKRKVTTNKC LEQVSQLQGL WRRFNRPLLK QQ.
Purity	> 97.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized TSLP in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Thymic Stroma...