


Product Datasheet

Human GRO-alpha protein (orb80128)

Description	GRO Alpha Human Recombinant produced in E.coli is a single, non-glycosylated, polypeptide chain containing...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	GRO a
Preservatives	The protein was lyophilized with no additives.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Ser-Val-Ala-Thr.
Purity	> 95.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized CXCL1 in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human GRO-alpha protein.