

Product Datasheet

Mouse MIP2 protein (orb80132)

Description	GRO-Beta Mouse Recombinant also caled mouse MIP-2 produced in E.coli is a single,non-glycosylated, p...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	mGRO b
Preservatives	The protein was lyophilized with no additives.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Protein quantitation was carried out by two independent methods1. UV spectroscopy at 280 nm using the absorbency value of 0.015 as the extinction coefficient for 0.1% (1mg/ml) solution. This value is calculated by the PC GENE computer analysis program of protein sequences (IntelliGenetics). 2. Analysis by RP-HPLC, using standard solution of GRO-b as Reference Standard.
Protein Sequence	The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Val-Val-Ala-Ser.
Purity	> 95.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized GRO-beta Mouse in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.