


Product Datasheet

Human PTH (7-34) protein (orb80169)

Description	Human PTH (7-34) protein
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	Human PTH (7-34)
Preservatives	The protein (1 mg/ml) was lyophilized after extensive dialyses against 20 mM PBS pH-7.0 and 4% mannitol.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	Leu-Met-His-Asn-Leu-Gly-Lys-His-Leu-Asn-Ser-Met-Glu-Arg-Val-Glu-Trp-Leu-Arg-Lys-Lys-Leu-Gln-Asp-Val-His-Asn-Phe.
Purity	> 97.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized PTH in sterile H ₂ O not less than 100 ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human PTH (7-34) pr...