

Product Datasheet

Bovine Enteropeptidase protein (orb80170)

Description	Enterokinase (rEK) Bovine Recombinant is the catalytic subunit of bovine enterokinase, which is expr...
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	bEnterokinase
Preservatives	Bovine EK in 50mM Tris-HCl, pH 8.0, 0.5M NaCl and 50% glycerol.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Application notes	One unit is defined as the amount of enzyme needed to cleave 50ug of fusion protein in 16 hours to 95% completion at 25C in buffer containing 25mM Tris-HCl, pH 7.6, 50mM NaCl, and 2mM CaCl2.
Expiration Date	6 months from date of receipt.