


Product Datasheet

Human MHC class chain-related gene A protein (orb80176)

Description	MICA Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing...
Conjugation	Unconjugated
Tested Applications	FA, HPLC, SDS-PAGE
Target	MICA
Preservatives	Lyophilized from concentrated (1mg/ml) solution containing no additives.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Protein Sequence	EPHSLRYNLTVLSWDGVSQSGFLAEVHLDGQPFLRYDRQKCRAPQGQWAED VLGNKTWDRETRDLTGNGKDLRMTLAHIKDQKEGLHSLQEIRVCEIHEDNSTRS SQHFYYDGELFLSQNLETEEWTPQSSRAQTLAMNVRNFLKEDAMKTKTHYHA MHADCLQELRRYLESGVLRRTVPPMVNTRSEASEGNITVTCRASSFYPRNIILT WRQDGVLSLHDTQQWGDVLPDGNQTYQTWWATRIRGEEQRFTCYMEHSGN HSTHPVPSGKVLVLQSH.
Purity	> 95.0% as determined by RP-HPLC and analysis by SDS-PAGE
Solubility (25°C)	It is recommended to reconstitute the lyophilized MICA in sterile H ₂ O not less than 100ug/ml, which can then be further diluted to other aqueous solutions.
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human MHC class cha...