

Product Datasheet

KDSA antibody (orb803594)

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Description | KDSA antibody |
| Species/Host | Rabbit |
| Reactivity | Bacteria |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant <i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2) <i>kdsA</i> protein |
| Target | <i>kdsA</i> |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q31KV0 |
| Expiration Date | 12 months from date of receipt. |