

Product Datasheet

MRAY antibody (orb803619)

Description	MRAY antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2) <i>mraY</i> protein
Target	<i>mraY</i>
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q31RY5
Expiration Date	12 months from date of receipt.