

Product Datasheet

RNA-directed RNA polymerase L antibody (orb805734)

Description	RNA-directed RNA polymerase L antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Drosophila melanogaster sigma virus (isolate Drosophila/USA/AP30/2005) (DMelSV) L protein
Target	L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C6S3N3
Expiration Date	12 months from date of receipt.