

Product Datasheet

Human SAA1 Protein (His Tag) (orb80598)

Description	Recombinant of Human SAA1 protein (His Tag)
Preservatives	The SAA1 solution contains 20mM Tris buffer(pH 8.0) and 10% glycerol.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MRSFFSFLGE AFDGARDMWR AYSDMREANY IGSDKYFHAR GNYDAAKRGP GGVWAAE AIS DARENIQRFF GHGAEDSLAD QAANEWGRSG KDPNHFRPAG LPEKY
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.