

## Product Datasheet

### Human TNFR2 Protein (His Tag) (orb80625)

<b>Description</b>	Recombinant of Human TNFR2 protein (His Tag)
<b>Preservatives</b>	TNFR2 protein is supplied in 20mM Tris HCl pH-8, 5mM EDTA and 50% glycerol.
<b>Form/Appearance</b>	Sterile Filtered clear solution.
<b>Storage</b>	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	Cytokines And Growth Factors
<b>Protein Sequence</b>	LPAQVAFTPYAPEPGSTCRLREYYDQTAQMCCSKCSPGQHAKVFCTKTSDTVCDSCEDSTY TQLWNWVPECLSCGSRCSDDQVETQACTREQNRICTCRPGWYCALSQEGCRLCAPLRK CRPGFGVARPGTETSDVVCKPCAPGTFSNTTSSDTCRPHQICNVVAIPGNASMDAVCTST SPT
<b>Purity</b>	Greater than 95.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	6 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)