

Product Datasheet

GE25944 antibody (orb807064)

| | |
|----------------------------|--|
| Description | GE25944 antibody |
| Species/Host | Rabbit |
| Reactivity | Drosophila |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Drosophila yakuba (Fruit fly) GE25944 protein |
| Target | GE25944 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | B4PUM2 |
| Expiration Date | 12 months from date of receipt. |