

Product Datasheet

Ganirelix Protein (orb80715)

Description	Ganirelix protein
Preservatives	The Ganirelix hormone (0.5mg/ml) contains 0.1mg acetic acid and 23.5mg mannitol pH-5.
Form/Appearance	Sterile Filtered colorless liquid formulation.
Storage	Stability: Ganirelix should be stored between 2°C- 8°C at all time. DO NOT FREEZE
Note	For research use only
Application notes	Hormones
Purity	Greater than 98.0% as determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.
Expiration Date	12 months from date of receipt.