

Product Datasheet

Human CSTA Protein (His Tag) (orb81122)

Description	Recombinant of Human CSTA protein (His Tag)
Preservatives	Cystatin-A (1mg/ml) in 20mM Tris-HCl buffer (pH8.0), 1mM DTT and 10% glycerol.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MIPGGLSEAK PATPEIQEIV DKVKPQLEEK TNETYGKLEA VQYKTQVVAG TNYIYKVRAG DNKYMHLKVF KSLPGQNE DL VLTGYQVDKN KDELDTGF
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.