

## Product Datasheet

### Human Myoglobin Protein (orb81226)

<b>Description</b>	Recombinant of Human Myoglobin protein
<b>Preservatives</b>	The sterile solution contains phosphate-buffered saline (pH 7.4) and 0.05% NaN <sub>3</sub> .
<b>Form/Appearance</b>	Sterile Filtered brownish solution.
<b>Storage</b>	Stability: Myoglobin should be stored at 4°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please do not freeze
<b>Note</b>	For research use only
<b>Application notes</b>	Recombinant & Natural Proteins
<b>Purity</b>	Greater than 95.0% as determined by SDS-PAGE.
<b>Source</b>	Escherichia Coli
<b>Expiration Date</b>	12 months from date of receipt.