

Product Datasheet

Human Myoglobin Protein (orb81227)

Description	Human Myoglobin protein
Preservatives	The protein solution is in 0.05M phosphate buffer pH 7.5 containing 0.15M NaCl and 0.09% NaN ₃ . Filtered through a 0.2µM membrane.
Form/Appearance	Sterile Filtered red solution.
Storage	Stability: Human Myoglobin should be stored at 2-8°C
Note	For research use only
Application notes	Human Virus Test: Starting material donor tested and certified negative for HIV I & II antibodies, Hepatitis B surface antigen, and Hepatitis C antibodies
Purity	Greater than 96.0%.
Source	Human Cardiac Tissues
Expiration Date	12 months from date of receipt.