

Product Datasheet

DBAT antibody (orb812351)

Description	DBAT antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant <i>Taxus cuspidata</i> (Japanese yew) 10-deacetylbaccatin III 10-O-acetyltransferase protein
Target	10-deacetylbaccatin III 10-O-acetyltransferase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9M6E2
Expiration Date	12 months from date of receipt.