

Product Datasheet

CYSK antibody (orb812465)

Description	CYSK antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Citrullus lanatus (Watermelon) (Citrullus vulgaris) Cysteine synthase protein
Target	Cysteine synthase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q43317
Expiration Date	12 months from date of receipt.