

Product Datasheet

Gapdhs antibody (orb815087)

Description	Gapdhs antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rattus norvegicus (Rat) Gapdhs protein
Target	Gapdhs
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9ESV6
Expiration Date	12 months from date of receipt.