

Product Datasheet

HIV-1 p24 Core Protein (orb81691)

Description	Recombinant of HIV-1 p24 Core protein
Preservatives	50mM Tris, pH 8.5 and 8M urea.
Form/Appearance	Sterile filtered colorless clear solution.
Storage	Stability: HIV-1 p24°Core although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles
Note	For research use only
Application notes	Viral Antigens- Human Immunodeficiency Virus
Purity	Greater than 95.0% as determined by HPLC analysis and SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.