

Product Datasheet

HIV-1 gag p17,p24, gp120 Protein (orb81718)

Description	Recombinant of HIV-1 gag p17,p24, gp120 protein
Preservatives	8M urea, 25mM Tris-HCl pH 8.0, 5mM b-mercaptoethanol.
Form/Appearance	Sterile filtered colorless clear solution.
Storage	Stability: HIV-1 gag p17,p24, gp120 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles
Note	For research use only
Application notes	Applications: HIV-1 gag p17,p24, gp120 antigen is suitable for ELISA and Western blots, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems
Purity	Greater than 95.0% as determined by HPLC analysis and SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.