

Product Datasheet

UL83 antibody (orb824341)

Description	UL83 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human cytomegalovirus (strain Towne) (HHV-5) (Human herpesvirus 5) UL83 protein
Target	UL83
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18139
Expiration Date	12 months from date of receipt.